

# W-210C Thermostatic Digital Laboratory Oil Bath



## Basic Info

- Model NO.: W-210C
- Structure: Desktop
- Certification: CE, ISO, SGS
- Type: Heating Type
- Trademark: JY
- Origin: Changzhou, China
  
- Customized: Customized
- Material: Stainless Steel
- Application: School, Hospital, Lab
- Size of Tank: 26X16cm
- Specification: 1PC in 1 carton
- HS Code: 8419899090

## Product Description

Thermostatic Digital Laboratory Oil Bath

1. LED, small size, RT-100d, hi-power
2. ISO13485/14001/9001 Certified
3. Since 1984

<b>Specifications of W-210C Thermostatic Digital Laboratory Oil Bath</b>	
Power Supply	AC 220V/110V 50Hz/60Hz
Temperature control range	Room temperature to 300 celsius
Fluctuate of Constant Temperature	± 1Celsius
Temperature coherence	1Celsius
Digital temperature control	Yes
Heating power	1000W
Working room	Φ 26× 16cm
Workarea: dia26cm*16cm, easy to use, rich function configuration, high quality, wide applicability. OEM service is offered. Match with a rotary evaporator.	