

Measure



WQS Vibrator With Sieve

Characteristic:

- * Particle analysis.
- * Widely used in the field of pharmaceutical industry, food, cosmetics, grain.
- * Advanced integrated circuit control
- * Easy to operate
- * The sieves is chosen by user (200,100,80,60,40meshsieve)

Specification:

- * Vibration direction: Vertical vibration
- * Measurement range: 325mesh
- * Vibration Frequency: 3000times/min,6000times/min
- * Amplitude selection: 0-3mm,continuousadjustment
- * Vibration mode: fine vibration; intermittent vibration; continuation vibration
- * Weight(Net):25kg