

X116(LED)



Specification		Model				
		X116(LED)	XSP-116HD	XSP-116HDS	XSP-116HB	XSP-116HBS
Viewing Head	Monocular Inclined at 38 °	•				
	Monocular Inclined at 45 °, 360° Rotatable		•	•	•	•
Eyepiece	WF10X	•	•	•	•	•
Objective	Achromatic 4X, 10X, 40X(s)	•	•	•	•	•
Stage	Plain stage with slide clips Stage size: 90mmX90mm	•				
	Plain stage with slide clips Stage size: 95mmX95mm		•	•	•	•
Condenser	0.65 Condenser with disc diaphragm	•	•	•	•	•
Focusing	Coarse Focusing adjustment	•	•		•	•
	Coaxial coarse & fine focusing adjustment			•		
Light Source	Transmitted Light LED, 110V/220V		•	•		
	Incident and transmitted Light LED, 110V/220V	•			•	•
Optional Accessory	Eye piece: WF16X, WF20X, P16X Objective: 20X, 60X(s)					