

XJL-03



- Small cubage, lightweight, easy operation and the wide range of accessories.
- With bright field, dark field and polarization observation system.
- Widely used for field metallurgical analysis and research in factories. Also used for image and video analysis.

Specification :

		XJL-03
Eyepiece	10×/16	●
	Eyepiece with scale of crosshair 10×	●
	12.5×/14	○
	16×/12	○
Plan Achromatic Objective	4×/0.10 WD 52.3mm	●
	10×/0.20 WD 17.3mm	●
	25×/0.40 WD 0.90mm	●
	40×/0.65 WD 0.40mm	●
	63×/0.85 WD 0.20mm	●
Zooming system	0.8×、1×、1.25×	●
Adapter for Dark Field	Adapter for objective 10×	●
	Adapter for objective 25×	●
Viewing Head	Sliding binocular viewing head	●
Kohler Incident Light	12V/50W Halogen light, Center, and brightness adjustable	●
	Polarizer and analyzer	●
Stage	Cycle-mechanical stage Φ160mm, Center adjustable, moving range 12mm×16mm	●
	小平台 : Φ10, Φ13, Φ15, Φ20, Φ21, Φ27 Stage: Φ10, Φ13, Φ15, Φ20, Φ21, Φ27	●
Accessories	Photo/Video attachment	○
	Micrometer-scale 0.01mm	●

Note: ● Standard Outfit, ○ Optional