

## XS-912



Achromatic objective: 4X,10X,40X(S),100X(S,OIL)

Eyepiece: WF 10X/18mm

Eyepiece head: Hinged Trinocular Head

Nosepiece: Quadruple

Stage: Double Layer Mechanical Stage Size 130X125mm

Focusing: Coaxial Coarse and Fine Adjustment, Focusing Range 30mm, Focusing Interval 0.002mm.

Condenser: Abbe NA=1.25 with Iris Diaphragm & Filter

Illumination: 6V/20W Halogen Lamp