

## **XSP-106**



Achromatic objective: 4X,10X,40X(S),100X(S,OIL)

Eyepiece: WF 10X

Eyepiece head: Hinged Binocular Head

Nosepiece: Quadruple

Stage: Double Layer Mechanical

Stage Size 125X115mm

Focusing: Coaxial Coarse and Fine Adjustment, Focusing Range 20mm, Focusing Interval 0.002mm.

Condenser: Abbe NA=1.25 with Iris Diaphragm & Filter

Illumination: 6V/20W Halogen Lamp