

XSP-116 Series



XSP-116



XSP-116L



XSP-116A



XSP-116AL



XSP-116D



XSP-116DL

Specification		Model					
		XSP-116	XSP-116L	XSP-116A	XSP-116AL	XSP-116D	XSP-116DL
Viewing Head	Monocular inclined at 30°	•	•	•	•	•	•
	360° Rotatable					•	•
Eyepiece	WF10X	•	•	•	•	•	•
Objective	Achromatic 4X、10X、40X (s)	•	•	•	•	•	•
Stage	Plain stage with slide clips	•	•	•	•	•	•
	Stage size: 90mmX90mm	•	•	•	•	•	•
Condenser	N.A.0.65 Condenser with disc diaphragm	•		•		•	•
Focusing	Coarse Focusing Adjustment		•		•		•
Light Source	Reflective Mirror	•		•		•	
Optional Accessory	Eyepiece: WF16X、P16X Objective: 20X、60X (s) LED illumination						