

XSP-116H



XSP-116H



XSP-116HD



XSP-116HB



XSP-116HBS



XSP-116HDS

Specification		Model				
		XSP-116H	XSP-116HD	XSP-116HDS	XSP-116HB	XSP-116HBS
Viewing Head	Monocular Inclined at 38 °	•				
	Monocular Inclined at 45 °, 360°Rotatable		•	•	•	•
Eyepiece	WF10X	•	•	•	•	•
Objective	Achromatic 4×、10×、40×(s)	•	•	•	•	•
Stage	Plain stage with slide clips Stage size: 90mmX90mm	•				
	Plain stage with slide clips Stage size: 95mmX95mm		•	•	•	•
Condenser	0.65 Condenser with disc diaphragm	•	•	•	•	•
Focusing	Coarse Focusing adjustment	•	•		•	•
	Coaxial coarse & fine focusing adjustment			•		
Light Source	Transmitted Light LED, 110V/220V		•	•		
	Incident and transmitted Light LED, 110V/220V	•			•	•
Optional Accessory	Eyepiece: WF16×、WF20×、P16× Objective: 20×、60×(s)					