

XSP-116NL Series



XSP-116NL



XSP-116B



XSP-116F



XSP-116FB

Specification		Model			
		XSP-116NL	XSP-116B	XSP-116F	XSP-116FB
Viewing Head	Monocular inclined at 45°	•	•	•	•
	360° Rotatable		•	•	•
Eyepiece	WF10X	•	•	•	•
Objective	Achromatic 4X、10X、40X (s)	•	•	•	•
Stage	Plain stage with slide clips Stage size: 90mmX90mm	•	•	•	•
Condenser	Single lens with Disc diaphragm	•	•	•	•
Focusing	Coarse Focusing Adjustment	•		•	•
	Coaxial coarse fine focusing adjustment Fine focusing scale value 0.002mm		•		•
Light Source	Incandescent lamp 220V/20W or 110V/20W	•	•	•	•
Optional Accessory	Eyepiece: WF16X、WF20X Objective: 20X、60X (s) LED illumination				