

XSP-117 Series



XSP-117



XSP-117F



XSP-117-100



XSP-117-200

Specification		Model			
		XSP-117	XSP-117F	XSP-117-100	XSP-117-200
Viewing Head	Monocular Inclined at 45°, 360° Rotatable	●	●	●	
	Sliding Binocular Head, Inclined at 45° (55mm-75mm)				●
Eyepiece	WF10×/18mm	●	●	●	●
Objective	Achromatic 4×, 10×, 40×(s)	●	●		
	Achromatic 4×, 10×, 40×(s), 100×(s) oil			●	●
Stage	Plain stage with slide clips Stage size: 120mm×120mm	●	●		
	Mechanical stage Moving range: 60mm×30mm			●	●
Condenser	N.A.0.65 Condenser with disc diaphragm	●			
	N.A.0.65 Condenser with iris diaphragm		●		
	N.A.1.25 Abbe Condenser with iris diaphragm and filter			●	●
Focusing	Separated coarse & fine focusing adjustment	●	●		
Light Source	Incandescent lamp 220V/20W or 110V/20W Or Fluorescent lamp 220V/5W or 110V/5W	●	●		
	Halogen bulb 6V/20W, AC220V/110V			●	●
Optional Accessory	Eyepiece: WF16×, P16× Objective: 20×, 60×(s) LED illumination				