

XSP-128 Series



XSP-128A



XSP-128B



XSP-128C



XSP-128T

Specification		Model				
		XSP-128	XSP-128A	XSP-128B	XSP-128C	XSP-128T
Viewing Head	Monocular Inclined at45°, 360°Rotatable	•	•	•		
	Compensation Free Binocular Head, Inclined t30°(55mm-75mm)				•	
	Compensation Free Trinocular Head, Inclined at30°(55mm-75mm)					•
Eyepiece	WF10×/18mm	•	•	•	•	•
Objective	Achromatic 4×、 10×、 40×(s)	•	•			
	Achromatic4×、 10×、 40×(s)、 100×(s)oil			•	•	•
Stage	Plain stage with slide clips Stage size:120mm×125mm	•	•			
	Mechanical stage Moving range:60mm×30mm			•	•	•
Condenser	0.65 Condenser with disc diaphragm	•				
	0.65 Condenser with iris diaphragm		•			
	N.A.1.25 Abbe Condenser with iris diaphragm And filter			•	•	•
Focusing	Separated coarse & fine focusing adjustment	•	•	•		•
Light Source	Incandescent lamp 220V/20W or 110V/20W	•				
	Fluorescent lamp 220V/5W or 110V/5W		•			
	Halogen bulb 6V/20W,AC 220V/110V Adjustable light intensity			•	•	•
Optional Accessory	Eyepiece:WF16×、 WF20×、 Pl6× Objective:20×、 60×(s) Reflective Mirror, Demonstration(Teaching)Head Compensation Free Binocular & Trinocular Head(48mm-75mm) LED illumination					