

XSP-13 Series



Specification		Model		
		XSP-13	XSP-17	XSP-17B
Viewing Head	Vertical Tube	●	●	●
Eyepiece	Huygens H10X、H16X	●		
	Huygens H5X 、H10X、H16X		●	
	Huygens H6X 、H10X、H15X			●
Objective	Achromatic 4X、10X、40X (s)	●		
	Achromatic 4X、10X、40X (s) oil		●	
	Achromatic 5X、10X、45X			●
Stage	Round stage with oris diaphragm			●
	Plain stage with disc diaphragm	●		
	Moving stage measuring tape Moving range: 60mmX30mm		●	
Condenser	N.A.1.25 Abbe Condenser with iris diaphragm and filter		●	
Focusing	Separated Coarse & Fine Focusing Adjustment	●	●	●
Light Source	Reflective Mirror	●	●	●
Optional Accessory	Eyepiece H12.5X Electrical light source 220V15W or 110V/15W			