

XSP-52 Series



Specification		Model		
		XSP-52	XSP-52F	XSP-52FL
Viewing Head	Monocular Head Vertical Tube	•	•	•
Eyepiece	WF10X with the pointer	•	•	•
Objective	Achromatic 4X、10X、40X (s)	•	•	•
Stage	Plain stage with slide clips	•	•	•
Condenser	Single lens with Disc diaphragm	•	•	•
Focusing	Coarse Focusing Adjustment	•		
	Separated coarse & fine focusing adjustment		•	•
Light Source	Plan-concave Reflective Mirror	•	•	
	Electrical light source 220V15W or 110V/15W			•
Optional Accessory	Eyepiece: WF16X			