

## XSP-53 Series



Specification		Model			
		XTD-223	XTD-223A	XTD-223T	XTD-223AT
Viewing Head	Binocular head (55mm-75mm)	•		•	
	Trinocular head (55mm-75mm)		•		•
Eyepiece	Plan WF10 X/ 22mm	•	•	•	•
Objective	Zoom range:0.8X-5.6X Zoom ratio: 1:7	•	•	•	•
	Working distance:110mm				
Light	Incident illumination LED/3W Brightness adjustable	•	•	•	•
	Transmitted illumination LED/1W Brightness adjustable	•	•	•	•
Optional Accessory	Eyepiece: WF10/23mm, WF15/17mm, WF20/12.5mm Right light Auxiliary Objective: 0.5X、 2X Photography attachment and CCD C-mount 1.3Mega、 2.0 Mega、 3.0Mega、 5.0Mega、 10.0Mega CMOS digital camera eyepiece Color CCD Camera,700 TV lines				