

XSP-68L Series



XSP-68L



XSP-68F



XSP-68H



XSP-68T

| Specification | | Model | | | |
|--------------------|---|---------|---------|---------|---------|
| | | XSP-68L | XSP-68F | XSP-68H | XSP-68T |
| Viewing Head | Monocular Inclined at 45°, 360° Rotatable | • | • | • | |
| | Sliding Binocular Head, Inclined at 45° (55mm-75mm) | | | | • |
| Eyepiece | WF10×/18mm | • | • | • | • |
| Objective | Achromatic 4×、10×、40×(s) | • | • | | |
| | Achromatic 4×、10×、40×(s)、100×(s)oil | | | • | • |
| Stage | Plain stage with slide clips Stage size: 120mm×110mm | • | • | | |
| | Mechanical stage Moving range: 60mm×30mm | | | • | • |
| Condenser | N.A.0.65 Condenser with disc diaphragm | • | | | |
| | N.A.0.65 Condenser with iris diaphragm | | • | | |
| | N.A.1.25 Abbe Condenser with iris diaphragm And filter | | | • | • |
| Focusing | Separated coarse & fine focusing adjustment | • | • | • | • |
| Light Source | Incandescent lamp 220V/20W or 110V/20W | • | • | • | • |
| Optional Accessory | Eyepiece: WF16×、WF20×、P16× Objective: 20×、60×(s) Demonstration (Teaching) Head WP-T or WP-V LED illumination | | | | |