

XSZ-120 Series



XSZ-120



XSZ-120E

Specification		Model	
		XSZ-120	XSZ-120E
Viewing Head	Compensation Free Binocular Head, Inclined at 30° (55mm-75mm)	•	
	Compensation Free Trinocular Head, Inclined at 30°(55mm-75mm)		•
Eyepiece	WF10×/18mm	•	•
Objective	Achromatic 4×、 10×、 40×(s)、 100×(s)Oil	•	•
Stage	Double layer mechanical stage		
	Stage Size:160mm×140mm	•	•
	Moving Range:80mm-50mm		
Condenser	N.A.1.25 Abbe condenser with iris diaphragm and filter	•	•
Focusing	Coaxial coarse & fine focusing adjustment Fine focusing scale value 0.002mm	•	•
Light Source	Halogen bulb 6V20W or 6V15W, AC 220V/110V Adjustable light source	•	•
Optional Accessory	Eyepiece: WF16×、 WF20×、 P16× Objective: 20×、 60(s) Semi-plan Objective:4×、 10×、 20×、 40×(s)、 60×(s)、 100×(s)oil Plan Objective:4×、 10×、 20×、 40×(s)、 60×(s)、 100×(s)oil Reflective Mirror, Demonstration(Teaching) Head Darkfield and Phase Contrast device 1.3Mega、 2.0Mega、 3.0Mega、 5.0Mega pixels CMOS Digital camera eyepiece Photography attachment and CCD Adapter 0.5×、 0.57×、 0.75× Compensation Free Binocular & Trinocular Head (48mm-75mm) LED illumination		