

## XSZ-127 Series



**XSZ-127**



**XSZ-127B**



**XSZ-127BT**

Specification		Model		
		XSZ-127	XSZ-127B	XAZ-127BT
Viewing Head	Monocular Inclined at 45°, 360° Rotatable	●		
	Compensation Free Binocular Head, Inclined at 30° (55mm-75mm)		●	
	Compensation Free Trinocular Head, Inclined at 30° (55mm-75mm)			●
Eyepiece	WF10×/18mm	●	●	●
Objective	Achromatic 4×、10×、40×(s)、100×(s)Oil	●	●	●
Stage	Double layer mechanical stage			
	Stage Size: 140mm×140mm	●	●	●
	Moving Range: 75mm×50mm			
Condenser	N.A.1.25 Abbe condenser with iris diaphragm and filter	●	●	●
Focusing	Coaxial coarse & fine focusing adjustment Fine focusing scale value 0.002mm	●	●	●
Light Source	Halogen bulb 6V20W or 6V15W, AC 220V/110V Adjustable light intensity	●	●	●
Base Frame	Black decorating panel	●	●	●
Optional Accessory	Eyepiece: WF16×、WF20×、P16× Objective: 20×、60(s) Semi-plan Objective: 4×、10×、20×、40×(s)、60×(s)、100×(s)oil Plan Objective: 4×、10×、20×、40×(s)、60×(s)、100×(s)oil Reflective Mirror, Demonstration(Teaching) Head Darkfield and Phase contrast device 1.3Mega、2.0Mega、3.0Mega、5.0Mega pixels CMOS Digital camera eyepiece Photography attachment and CCD Adapter 0.5×、0.57×、0.75× Compensation Free Binocular & Trinocular Head (48mm-75mm) LED illumination			