

XSZ-127DA Series



XSZ-127DA



XSZ-127E



XSZ-127D



XSZ-127DE

Specification		Model			
		XSZ-127DA	XSZ-127E	XSZ-127D	XSZ-127DE
Viewing Head	Sliding Binocular Head, Inclined at 45 (55mm-75mm)	•		•	
	Sliding Trinocular Head, Inclined at 45 (55mm-75mm)		•		•
Eyepiece	WF10×/18mm	•	•	•	•
Objective	Achromatic 4×、10×、40×(s)、100×(s) Oil	•	•	•	
Stage	Double layer mechanical stage				
	Stage Size:140mm140mm	•	•	•	•
	Moving Range:75mm×50mm				
Condenser	N.A.1.25 Abbe condenser with iris diaphragm and filter	•	•	•	•
Focusing	Coaxial coarse & fine focusing adjustment in focusing scale value 0.002mm	•	•	•	•
Light Source	Halogen bulb 6V20W or 6V15W, AC 220V/110V adjustable light intensity	•	•	•	•
Base Frame	Black decorating panel			•	•
Optional Accessory	Eyepiece: WF16×、WF20×、P16× Objective: 20×、60(s)、 Semi-plan Objective:4×、10×、20×、40×(s)、60×(s)、100×(s)oil Plan Objective:4×、10×、20×、40×(s)、60×(s)、100×(s)oil Reflective Mirror, Demonstration(Teaching) Head Darkfield, and Phase contrast device 1.3Mega、2.0Mega、3.0Mega、 5.0Mega pixels CMOS Digital camera eyepiece Photography attachment and CCD Adapter 0.5×、0.57×、0.75× LED illumination				