

## XSZ-135A Series



Specification		Model		
		XSZ-135A	XSZ-135B	XAZ-135T
Viewing Head	Monocular Inclined at 45°, 360° Rotatable	•		
	Compensation Free Binocular Head, Inclined at 30° (55mm-75mm)		•	
	Compensation Free Trinocular Head, Inclined at 30° (55mm-75mm)			•
Eyepiece	WF10x/18mm	•	•	•
Objective	Achromatic 4x, 10x, 40x(s), 100x(s) Oil	•	•	•
Stage	Double layer mechanical stage			
	Stage Size: 140mm x 140mm	•	•	•
	Moving Range: 70mm x 50mm			
Condenser	N.A. 1.25 Abbe condenser with iris diaphragm and filter	•	•	•
Focusing	Coaxial coarse & fine focusing adjustment Fine focusing scale value 0.002mm Focusing Range 30mm	•	•	•
Light Source	Halogen bulb 6V20W or 6V15W, AC 220V/110V Adjustable light source with Kohler illumination	•	•	•
Optional Accessory	Eyepiece : WF16x, WF20x, P16x Objective : 20x, 60(s), Semi-plan Objective: 4x, 10x, 20x, 40x(s), 60x(s), 100x(s) oil Plan Objective: 4x, 10x, 20x, 40x(s), 60x(s), 100x(s) oil, Reflective Mirror, Demonstration (Teaching) Head Darkfield and Phase contrast device 1.3Mega, 2.0Mega, 3.0Mega, 5.0Mega pixels CMOS Digital camera eyepiece Photography attachment and CCD Adapter 0.5x, 0.57x, 0.75x Compensation Free Binocular & Trinocular Head (48mm-75mm) LED illumination			