

XSZ-136A Series



XSZ-136A



XSZ-136B



XSZ-136E

Specification		Model		
		XSZ-136A	XSZ-136B	XSZ-136E
Viewing Head	Monocular Inclined at 45°, 360° Rotatable	●		
	Compensation Free Binocular Head, Inclined at 30° (55mm-75mm)		●	
	Compensation Free Trinocular Head, Inclined at 30° (55mm-75mm)			●
Eyepiece	WF10×/18mm	●	●	●
Objective	Achromatic 4×, 10×, 40×(s), 100×(s) oil	●	●	●
Stage	Double layers mechanical stage			
	Stage size: 135mm×135mm	●	●	●
	Moving Range: 70mm×30mm			
Condenser	N.A.1.25 Abbe Condenser with iris diaphragm and filter	●	●	●
Focusing	Coaxial coarse & fine focusing adjustment Fine focusing scale value 0.002mm	●	●	●
Light Source	Halogen bulb 6V/20W, AC 220V/110V	●	●	●
Optional Accessory	Eyepiece: WF16×, WF20× Objective: 20×, 60×(s) Compensation Free Binocular & Trinocular Head (48mm-75mm) LED illumination			