

XSZ-137 Series



Specification		Model		
		XSZ-137	XSZ-137E	XSZ-137F2
Viewing Head	Compensation Free Binocular Head, Inclined at 30° (55mm-75mm)	•		
	Compensation Free Trinocular Head, Inclined at 30° (55mm-75mm)		•	
	Demonstration(Teaching) binocular Head			•
Eyepiece	WF10x/18mm	•	•	•
Objective	Achromatic 4x、 10x、 40x(s)、 100x(s)Oil	•	•	•
Stage	Double layer mechanical stage			
	Stage Size:140mm×140mm	•	•	•
	Moving Range:70mm×50mm			
Condenser	N.A.1.25 Abbe condenser with iris diaphragm and filter	•	•	•
Focusing	Coaxial coarse & fine focusing adjustment Fine focusing scale value 0.002mm	•	•	•
Light Source	Halogen bulb 6V20W or 6V15W, AC 220V/110V	•	•	•
	Adjustable light source with Kohler illumination			
Optional Accessory	Eye piece : WF16x、 WF20x、 P16x Objective : 20x、 60(s) Semi-plan Objective:4x、 10x、 20x、 40x(s)、 60x(s)、 100x(s)oil Plan Objective:4x、 10x、 20x、 40x(s)、 60x(s)、 100x(s)oil、 Reflective Mirror, Demonstration(Teaching) Head Dark field and Phase contrast device 1.3Mega、 2.0Mega、 3.0Mega、 5.0Mega pixels CMOS Digital camera eyepiece Photography attachment and CCD Adapter 0.5x、 0.57x、 0.75x Compensation Free Binocular & Trinocular Head (48mm-75mm) LED illumination			