

XSZ-138 Series



XSZ-138



XSZ-138E

Specification		Model	
		XSZ-138	XSZ-138E
Viewing Head	Compensation Free Binocular Head, Inclined at 30°(55mm-75mm)	●	
	Compensation Free Trinocular Head, Inclined at 30°(55mm-75mm)		●
Eyepiece	WF10×/18mm	●	●
	WF10×/22mm		
Objective	Infinity Achromatic 4×、10×、40×(s)、100×(s)Oil		
	Achromatic 4×、10×、40×(s)、100×(s)Oil	●	●
Stage	Double layer mechanical stage		
	Stage Size:150mm×150mm	●	●
	Moving Range:80mm×50mm		
Condenser	N.A.1.25 Abbe condenser with iris diaphragm and filter	●	●
Focusing	Coaxial coarse & fine focusing adjustment Fine focusing scale value 0.002mm Focusing Range 30mm	●	●
Light Source	Halogen bulb 6V20W or 6V15W, AC 220V/110V Adjustable light source with Kohler illumination	●	●
Optional Accessory	Eyepiece: WF16×、WF20×、P16× Objective: 20×、60(s) Semi-plan Objective:4×、10×、20×、40×(s)、60×(s)、100×(s)oil Plan Objective:4×、10×、20×、40×(s)、60×(s)、100×(s)oil Reflective Mirror, Demonstration(Teaching) Head Dark field and Phase contrast device 1.3Mega、2.0Mega、3.0Mega、5.0Mega pixels CMOS Digital camera eyepiece Photography attachment and CCD Adapter 0.5×、0.57×、0.75× Compensation Free Binocular & Trinocular Head (48mm-75mm) LED illumination		