

XSZ-139 Series



XSZ-139



XSZ-139T



XSZ-139A



XSZ-139AT

Specification		Model			
		XSZ-139	XSZ-139T	XSZ-139A	XAZ-139AT
Viewing Head	Compensation Free Binocular Head, Inclined at 30° (55mm-75mm)	•		•	
	Compensation Free Trinocular Head, Inclined at 30° (55mm-75mm)		•		•
Eyepiece	WF10×/18mm	•	•		
	WF10×/22mm			•	•
Objective	Infinity Achromatic 4×、10×、40×(s)、100×(s)Oil			•	•
	Achromatic 4×、10×、40×(s)、100×(s)Oil	•	•		
Stage	Double layer mechanical stage				
	Stage Size:150mm×150mm	•	•	•	•
	Moving Range:80mm×50mm				
Condenser	N.A.1.25 Abbe condenser with iris diaphragm and filter	•	•	•	•
Focusing	Coaxial coarse & fine focusing adjustment Fine focusing scale value 0.002mm Focusing Range 30mm	•	•	•	•
Light Source	Halogen bulb 6V20W or 6V15W, AC 220V/110V Adjustable light source with Kohler illumination	•	•	•	•
Optional Accessory	Eyepiece×、WF20×、P16× Objective: 20×、60(s)、 Semi-plan Objective:4×、10×、20×、40×(s)、60×(s)、100×(s)oil Plan Objective:4×、10×、20×、40×(s)、60×(s)、100×(s)oil Reflective Mirror, Demonstration(Teaching) Head Dark field and Phase contrast device 1.3Mega、2.0Mega、3.0Mega、5.0Mega pixels CMOS Digital camera eyepiece Photography attachment and CCD Adapter 0.5×、0.57×、0.75× Compensation Free Binocular & Trinocular Head (48mm-75mm) LED illumination				