

XSZ-152A Series



Features:

The series microscope is widely used in Biology, Histology, Pathology, Medicament, chemistry research and Clinical examination. It can also be used for education in Colleges and Schools.

Specification		Model		
		XSZ-152A	XSZ-152B	XSZ-152T
Viewing head	Monocular Head. Inclined at 45°	●		
	Compensation Free Binocular Head .Inclined at 30°		●	
	Compensation Free Trinocular Head Inclined at 30°			●
Interpapillary Distance	55mm -75mm		●	●
Eyepiece	WF10X	●	●	●
Objective	Achromatic objectives:4X,10X,40X(S),100X(S) Oil	●	●	●
Stage	Double layers mechanical stage	●	●	●
	Stage size:125mmX 130mm	●	●	●
	Moving range: 75mm × 50mm	●	●	●
Condenser	Swing type N. A.1.25 Abbe condenser with iris diaphragm & filter	●	●	●
Focusing	Coaxial coarse & fine focusing adjustment with rack and pinion mechanism Fine focusing scale Value 0.002mm	●	●	●
Collector	High brightness light collector	●	●	●
Light source	Halogen bulb 6V/20W. AC85-230V	●	●	●
	Adjustable brightness			
Optional Accessory	Eyepieces: WF16x,WF20x,P16x Plan-Objectives: 4X,10X,20X(S),40X(S),60X(S),100X(S)Oil Semi-plan objectives: 4X,10X,20X(S),40X(S),60X(S),100X(S)Oil Photo and CCD attachment			