

## XSZ-157 Series



**XSZ-157B**

**XSZ-157E**

Specification		Model		
		XSZ-157	XSZ-157B	XAX-157E
Viewing Head	Monocular Inclined at 45°, 360° Rotatable	•		
	Sliding Binocular Head, Inclined at 45° (55mm-75mm)		•	
	Sliding Trinocular Head, Inclined at 45° (55mm-75mm)			•
Eyepiece	WF10×/18mm	•	•	•
Objective	Achromatic 4×, 10×, 40×(s), 100×(s) Oil	•	•	•
Nosepiece	Quadruple Nosepiece	•	•	•
Stage	Double layer mechanical stage			
	Stage Size: 140mm×132mm	•	•	•
	Moving Range: 70mm×30mm			
Condenser	N.A. 1.25 Abbe condenser with iris diaphragm and filter	•	•	•
Focusing	Coaxial coarse & fine focusing adjustment Fine focusing scale value 0.002mm	•	•	•
Light Source	Halogen bulb 6V20W or 6V15W, AC 220V/110V Adjustable light intensity	•	•	•
Optional Accessory	Eyepiece : WF16×, WF20×, P16× Objective : 20×, 60(s), Semi-plan Objective: 4×, 10×, 20×, 40×(s), 60×(s), 100×(s) oil Plan Objective: 4×, 10×, 20×, 40×(s), 60×(s), 100×(s) oil, Reflective Mirror, Demonstration (Teaching) Head Dark field and Phase contrast device 1.3Mega, 2.0Mega, 3.0Mega, 5.0Mega pixels CMOS Digital camera eyepiece Photography attachment and CCD Adapter 0.5×, 0.57×, 0.75× Compensation Free Binocular & Trinocular Head (48mm-75mm) LED illumination Quintuple Nose-piece			