

## XSZ-157A Series



**XSZ-157H**



**XSZ-157T**

Specification		Model		
		XSZ-157A	XSZ-157H	XAX-157T
Viewing Head	Monocular Inclined at 45°, 360° Rotatable	•		
	Compensation Free Binocular Head, Inclined at 30° (55mm-75mm)		•	
	Compensation Free Trinocular Head, Inclined at 30° (55mm-75mm)			•
Eyepiece	WF10×/18mm	•	•	•
Objective	Achromatic 4×, 10×, 40×(s), 100×(s) Oil	•	•	•
Nosepiece	Backward Quadruple Nosepiece	•	•	•
Stage	Double layer mechanical stage			
	Stage Size: 140mm×132mm	•	•	•
	Moving Range: 70mm×30mm			
Condenser	N.A.1.25 Abbe condenser with iris diaphragm and filter	•	•	•
Focusing	Coaxial coarse & fine focusing adjustment Fine focusing scale value 0.002mm	•	•	•
Light Source	Halogen bulb 6V20W or 6V15W, AC 220V/110V	•	•	•
	Adjustable light intensity			
Optional Accessory	Eyepiece : WF16×, WF20×, P16× Objective : 20×, 60(s) Semi-plan Objective: 4×, 10×, 20×, 40×(s), 60×(s), 100×(s) oil Plan Objective: 4×, 10×, 20×, 40×(s), 60×(s), 100×(s) oil Reflective Mirror, Demonstration (Teaching) Head Dark field and Phase contrast device 1.3Mega, 2.0Mega, 3.0Mega, 5.0Mega pixels CMOS Digital camera eyepiece Photography attachment and CCD Adapter 0.5×, 0.57×, 0.75× Compensation Free Binocular & Trinocular Head (48mm-75mm) LED illumination			