

## XSZ-158



### Features:

XSZ-158 microscope becomes a new landmark since Ningbo teaching instrument Co., Ltd was founded. It is diversification brings you enter into an epoch not only the performance but also the application. The module design can meet your needs on all normal experiment and scientific research. Its humanized configuration design, reasonable control device, and layout of every operation knob can keep your eyes, neck, shoulders, waist, arms, and wrist against tiredness. You can be released from heavy word due to the comfortable operation.

Specification		Standard outfit	Optional parts
Viewing head	Compensation Free Binocular Head, Inclined adjustment 0 ° -40 °	●	
	Compensation Free Trinocular Head ,Inclined adjustment 0 ° -40 °		○
Interpapillary Distance	55mm -75mm	●	
Eyepiece	WF10X/ 22mm	●	
	WF16X,WF20X,WF25X,P16X		○
Nosepiece	Quintuple nosepiece	●	
	Sextuple nosepiece		○
Objective	Infinity plan objectives:4X,10X,40X(S),100X(S) Oil	●	
	Infinity plan objectives:20X(S),60X(S)		○
Stage	Double layers mechanical stage	●	
	Stage size:180mmX 160mm		
	Moving range :75mmX 50mm		
Condenser	Swing type N. A.1.25 Abbe condenser with iris diaphragm & filter	●	
Focusing	Coaxial coarse & fine focusing adjustment with rack and pinion mechanism Fine focusing scale Value 0.002mm	●	
Collector	Kohler illumination	●	
Light source	Halogen bulb12V/30W 12V/50W or 12V/100W.AC85-230V	●	
	Adjustable brightness		
Optional Accessory	Fluorescent attachment		○
	Phase contrast attachment		○
	Dry & Immersion dark field condenser		○
	Simple polarization attachment		○
	3.0or1.3Mega pixel, CMOS electronic eyepiece		○