

XSZ-170



XSZ-170B



XSZ-170BT

Specification		Model			
		XSZ-170	XSZ-170T	XSZ-170B	XSZ-170BT
ViewingHead	Compensation Free Binocular Head, Inclined at 30° (48mm -75mm)	•		•	
	Compensation Free Binocular Head, Inclined at 30° (48mm -75mm)		•		•
Eyepiece	WF10 X 18mm	•	•		
	WF10 X 22mm			•	•
Objective	Achromatic Objective : 4X , 10X , 40X(S), 100X (S) Oil	•	•		
	infinity Achromatic Objective : 4X , 10X , 40X(S) , 100X (S) Oil			•	•
Stage	Double layer mechanical stage				
	Stage Size : 132mmx 140mm	•	•	•	•
	Moving Range : 70mmx 50mm				
Condenser	N. A.1.25 Abbe condenser with iris diaphragm and filter	•	•	•	•
Focusing	Coaxial coarse & fine focusing adjustment, fine focusing scale value 0.002mm focusing range 30mm	•	•	•	•
Light Source	LED, brightness adjustment	•	•	•	•
Optional Accessory	Eyepiece: WF16x,WF20x,P16x Achromatic Objective : 20X(S)、 60X(S) Plan objective : 4X,10X,20X(S),40X(S),60X(S),100X(S)Oil infinity plan objective: 4X,10X,20X(S),40X(S),60X(S),100X(S)Oil Dual viewing head Dark field and phase contrast attachment 1.3Mega、 3.0 Mega、 5.0Mega、 10.0Mega COMS digital camera eyepiece Photography attachment and CCD C-mount:0.5X、 0.57X、 0.75X CCD Camera,colur1/3 High resolution 700 TV lines				