

## XSZ-178M Series



Specification		Model		
		XSZ-178M	XSZ-178MB	XSZ-178MT
Viewing Head	Monocular head, inclined at 45 °, 360° rotatable	●		
	Compensation free binocular head, inclined at 30°(48mm-75mm)		●	
	Compensation free binocular head, inclined at 30°(48mm- 75mm)			●
Eyepiece	WF10×/18mm	●	●	●
Objective	Achromatic objective:4×、 10×、 40×(s)、 100×(s)Oil	●	●	●
Stage	Double layer mechanical stage			
	Stage size : 125mm×130mm	●	●	●
	Moving Range : 70mm×50mm			
Condenser	N.A.1.25 Abbe condenser with iris diaphragm & filter holder	●	●	●
Focusing	Coaxial coarse & fine focusing adjustment, fine focusing scale value 0.002mm focusing range 30mm	●	●	●
Light Source	LED, Brightness adjustable	●	●	●
Optional Accessory	Eye piece : WF16×、 WF20×、 P16× Achromatic Objective : 20× (s) 、 60×(s) Semi-plan objective : 4×、 10 ×、 20×、 40×(s)、 60×(s)、 100×(s)oil Plan objective : 4×、 10×、 20×、 40×(s)、 60×(s)、 100×(s)oil Dual viewing head Darkfield and phase contrast attachment 1.3Mega、 3.0Mega、 5.0Mega、 10.0Mega COMS digital camera eyepiece Photography attachment and CCD C-mount:0.5X、 0.57X、 0.75X CCD Camera, colour 1/3 High resolution 700 TV lines			