

XSZ-179



XSZ-179H



XSZ-179D



XSZ-179E

Specification		Model		
		XSZ-179H	XSZ-179D	XSZ-179E
Viewing Head	Monocular head, inclined at 45 °, 360° rotatable	•		
	Compensation free binocular head, inclined at 30°(48mm- 75mm)		•	
	Compensation free binocular head, inclined at 30°(48mm- 75mm)			•
Eyepiece	WF10×/18mm	•	•	•
Objective	Achromatic objective:4×、10×、40×(s)、100×(s)Oil	•	•	•
Stage	Double layer mechanical stage			
	Stage size : 135mm×135mm	•	•	•
	Moving Range : 70mm×50mm			
Condenser	N.A.1.25 Abbe condenser with iris diaphragm & filter holder	•	•	•
Focusing	Coaxial coarse & fine focusing adjustment, fine focusing scale value 0.002mm focusing range 30mm	•	•	•
Light Source	LED, Brightness adjustable	•	•	•
Optional Accessory	Eye piece : WF16×、WF20×、P16× Objective : 20× (s) 、60×(s) Reflector, Dual viewing head			