

## XSZ-208A Series



Specification		Model		
		XSZ-208A	XSZ-208	XSZ-208E
Viewing Head	Monocular Inclined at 45°, 360° Rotatable	●		
	Compensation Free Binocular Head, Inclined at 30° (55mm-75mm)		●	
	Compensation Free Trinocular Head, Inclined at 30° (55mm-75mm)			●
Eyepiece	WF10×/18mm	●	●	●
Objective	Achromatic 4×, 10×, 40×(s), 100×(s) oil	●	●	●
Stage	Double layers mechanical stage	●	●	●
	Stage size: 140mm×140mm			
	Moving Range: 75mm×50mm			
Condenser	N.A.1.25 Abbe Condenser with iris diaphragm and Filter	●	●	●
Focusing	Coaxial coarse & fine focusing adjustment Fine focusing scale value 0.002mm	●	●	●
Light Source	Halogen bulb 6V20W or 6V15W, AC220V/110V Adjustable light intensity	●	●	●
Optional Accessory	Eyepiece: WF16×, WF20×, P16× Objective: 20×, 60×(s) Semi-plan Objective: 4×, 10×, 20×, 40×(s), 60×(s), 100×(s) oil Plan Objective: 4×, 10×, 20×, 40×(s), 60×(s), 100×(s) oil Reflective Mirror, Demonstration (Teaching) Dark field and Phase contrast device 1.3Mega, 2.0Mega, 3.0Mega, 5.0Mega pixels Photography attachment and CCD Adapter 0.5×, Compensation Free Binocular & Trinocular Head LED illumination			