

XSZ-2108



Eyepiece:	WF 10X/18mm
Head:	Binocular head.
Objectives:	Achromatic 4X,10X,40X,100X
Nosepiece:	Quadruple
Condenser:	Abbe condenser, A=1.25(when immersed in oil)
Focusing:	Coaxial Coarse and Fine Adjustment, Focusing Range 22mm,Focusing Interval 0.002mm.
Stage: Size:	140X140mm, moving rang 75X45mm
Illumination:	6V20W, Halogen lamp.