

## **XSZ-2108A**



Optical system: Infinity Optical system

Eyepiece: WF 10X/18mm

Head: Binocular head.

Objectives: Infinity Plan 4X,10X,40X,100X

Nosepiece: Quadruple

Condenser: Abbe condenser, A=1.25(when immersed in oil)

Focusing: Coaxial Coarse and Fine Adjustment, Focusing Range 22mm,Focusing Interval 0.002mm.

Stage: Size: 140X140mm, moving rang 75X45mm

Illumination: 6V20W, Halogen lamp.