

## **XSZ-910**



Achromatic objective: 4X,10X,40X(S),100X(S,OIL)

Eyepiece: WF 10X/18mm

Eyepiece head: Hinged Binocular Head

Nosepiece: Quadruple

Stage: Double Layer Mechanical Stage Size 130X125mm

Focusing: Coaxial Coarse and Fine

Adjustment,Focusing Range: 30mm,Focusing Interval 0.002mm.

Condenser: Abbe NA= 1.25 with Iris Diaphragm &Filter

Illumination: 6V/20W Halogen Lamp