

## XTD-203 Series



XTD-203(A, B, C)



XTD-205(A, B, C)



XTD-206(A, B, C)



XTD-210(A, B, C)

This series microscope is an implemental stereo microscope. Because of eximious contrast .deep depth of field and high resolution, it is widely used in Agriculture, Jewelry, electronic industry, Precision mechanism, etc and also suitable for viewing the 3D object in both natural and engineering fields. It can be a standard tool for education and scientific research too.

Model	Viewing Head	Eyepieces	Objective	Working distance	Illumination	Optional parts
XTD -203A	Binocular Head Inclined at 45° 360° rotatable	SWF10X/20	1X 、 2X	80mm	Incident illumination 12V/10W or 12V/15W	Fluorescent ring lamp – 1.3Mega pixels electronic eyepiece – Eyepiece: SWF15X/15 SWF20X/11 -LED LED illumination
XTD-203B			1X 、 3X			
XTD -203C			2X 、 4X			
XTD -205A			1X 、 2X			
XTD-205B			1X 、 3X			
XTD -205C			2X 、 4X			
XTD -206A			1X 、 2X	110mm	Transmission illumination 12V/10W or 12V/15W	
XTD-206B			1X 、 3X			
XTD -206C			2X 、 4X			
XTD -210A			1X 、 2X			
XTD-210B			1X 、 3X			
XTD -210C			2X 、 4X			