

XZB-605



Model XZB-605 Comparison Microscope makes you view two images (left & right) from two objects at the same time. Adjusting two images to connect or lap over, then differentiate the minuteness difference between them, for instance, identifying the bullet's mark, tool's mark, fingerprint, stamp, letters, money, presswork, etc. It is an ideal instrument for Police and other departments, which need object comparison.

There're several different viewing state available on comparison head;

- Full left image, full right image, a part of left & right images, left & right superposition image.
- The total magnification of comparison head: 1X
- The viewing field of comparison head is more than 20mm.
- 195 series objectives are applicable.
- Biological microscope objectives 4X , 10X , 20X , 40X , 50X , 60X , 100X are suitable for this microscope.
- Monocular, Binocular or Trinocular Comparison head are for option.