

YB-1A VACUUM CONSTANT TEMPERATURE DRYING CASE



YB-1A VACUUM CONSTANT TEMPERATURE DRYING CASE

Applicable standard
National standard (Chinese Codex) /2015

Features:

- Operate easily with slight shake buttons.
- High vacuum seal.
- It can achieve an equal temperature as heating.
- Display vacuum by an analog meter and display temperature by a digital meter.
- Observe samples by the pressure-proof glass at all hours.

Specifications:

- Vacuum range : (0~0.1) Mp
- Temperature range: (30~150)°C
- Temperature Stability: $\pm 2^{\circ}\text{C}$
- Usage temperature: (15~35)°C
- Power: 220V/ 50Hz/ 150W
- Dimension: (360×360×250)mm³