

YD-1 TABLET HARDNESS TESTER



YD-1 TABLET HARDNESS TESTER

Tablet hardness testers are instruments for detecting breaking hardness of tablet.

Applicable standards:

Corporate standard(Tablet hardness tester)Q/12XQ0186-2010

Features:

The tablet is loaded and pressed artificially.

It has high accuracy press sensor and hardness data is displayed on LED.

It tests tablet hardness continuously.

Test data is displayed and locked automatically, so can test be reset and cycled.

Automation: auto-test, auto-diagnose, auto-alarm.

Specifications:

Hardness range: (2~199.9)N

Hardness accuracy Max.: *($\pm 0.05\%$)

Moving distance of probe: 20mm

Power: 220V/ 50Hz/ 10W or 110V/ 60Hz/ 10W

Dimension: (320*200*100)mm³