

SEED TESTING INSTRUMENT

Automatic Seed Counter

Model: SLY-C

●Characteristic

- ▲ Microcomputer automatic control, have touch-key, complete automated operation.
- ▲ Count speed is adjustable, high accuracy, No noise metal shell, and elegant appearance.
- ▲ Both round and long seeds are suitable
- ▲ Automatically adjust the grain of clearance

●Technical specification

Count Range: Small and middle sample: 0.7~4 mm×0.7~12mm
Big sample: 3~10mm×3~12 mm
Count Accuracy: Big and middle sample: 2/1000 small sample: 4/1000
Count Speed: ≥1000 /3min (it's based on level 2, if you want fast you can turn to level 3-7)
Count Capacity: 1~99999
Power: 220V 50HZ
Continuous working time: ≥ 5Hr

